

Apparatus A 3-mL and 5-mL conical vial, screw-cap centrifuge tube, two 1-mL plastic (or glass) syringes, Pasteur pipet with 0.5- and 1.0-mL calibration marks, Claisen adapter, drying tube, ice-water bath, and apparatus for magnetic stirring and flameless heating.

Setting Up Dry the conical vials, Claisen adapter, centrifuge tube, drying tube, and calibrated Pasteur pipet in an oven at 110 °C for at least 30 min. Do *not* put any plastic connectors or rubber O-rings in the oven, as they may melt or soften. Using gloves or tongs, remove the glassware from the oven and let it cool, preferably in a desiccator. Lubricate the joints and assemble the apparatus shown in Figure 19.1 by adding a spinvane to the 5-mL conical vial and then fitting the Claisen adapter to the vial. Finally, place a rubber septum and the drying tube on the Claisen adapter.

Optional Measures If an oven is not available, it will be necessary to dry the apparatus with a microburner or a heat gun. Assemble the apparatus as described above. *Be sure that no one in the laboratory is working with diethyl ether*, and then dry the assembled apparatus. Be careful heating the conical vial, as this heavy-walled vessel can easily crack from thermal shock. Also, do not overheat any plastic parts of the apparatus. Allow the apparatus to cool to room temperature.

Verify that there are no flames in the laboratory before continuing. Transfer about 3–4 mL of anhydrous diethyl ether to a screw-cap centrifuge tube and cap the tube. Use this during the remainder of the experiment whenever anhydrous diethyl ether is required. Weigh 0.05 g of magnesium turnings that have been freshly crushed with a heavy spatula or the end of a clamp. Remove the rubber septum and transfer the turnings to the reaction vial. Add one *small* crystal of iodine to the vial and replace the rubber septum. Add 0.2 mL of anhydrous diethyl ether to the vial using a dry 1-mL syringe inserted through the rubber septum. Prepare a solution of 0.24 mL of bromobenzene in 0.5 mL of anhydrous diethyl ether in a dry 3-mL conical vial. Swirl the solution to achieve homogeneity.




Figure 19.1
Microscale apparatus for
preparing and reaction of
Grignard reagents.

Reaction Stir the contents of the conical vial while warming them gently. Draw the ethereal solution of bromobenzene into the plastic syringe and cap the conical vial containing the bromobenzene solution. Insert the syringe needle through the rubber septum, and add a 0.1-mL portion of this solution onto the magnesium turnings; stir the mixture. *Hold* the plunger of the syringe to control the amount of solution added. If small bubbles form at the surface of the magnesium turnings or if the mixture becomes slightly cloudy or chalky, the reaction has started.

Optional Measures If the reaction does not start spontaneously, consult your instructor. The best remedy at this point is to continue warming the vial and add an additional crystal of iodine to the mixture. You may also add a small sample of phenylmagnesium bromide from a stock solution if it is available.

Once the reaction has started, continue heating the reaction mixture gently so the solvent refluxes slightly. Using a second plastic syringe inserted through the rubber septum, add another 0.5-mL portion of anhydrous diethyl ether to the reaction mixture and continue heating and stirring until the solvent is again refluxing. Add the remainder of the bromobenzene-ether solution *dropwise* to the stirred reaction mixture at a rate that is just fast enough to maintain a gentle reflux. If the reaction becomes too vigorous, reduce the rate of adding the ethereal solution of aryl halide and, if necessary, discontinue heating the vial. If the spontaneous boiling of the mixture slows, increase the rate of addition slightly. If the rate of reflux still does not increase, increase heating and maintain gentle reflux during the remainder of addition. *It is important that reflux be maintained throughout the addition of the bromobenzene-ether solution.* The addition should take about 3–5 min. Upon completion of the addition, place about 0.5 mL of anhydrous diethyl ether in the vial that contained the bromobenzene solution. Draw this solution into a syringe and add it in one portion to the reaction mixture. Continue heating the mixture with stirring under gentle reflux for 15 min. If necessary, add anhydrous diethyl ether to maintain a level *no lower* than the 2-mL mark on the conical vial. At the end of the reaction, the solution normally has a tan to brown, chalky appearance, and most of the magnesium will have disappeared, although residual bits of metal usually remain. Discontinue heating and allow the mixture to cool to room temperature.

Use the Grignard reagent as soon as possible after preparing it, following one of the procedures given in Section 19.4.